creos™ mucogain

Patient Information Leaflet

en	1
da	7
de	13
el	19
es	25
et	31
fi	37
fr	43
hr	49
hu	55
it	61
lt	67
lv	73
nl	79
no	85
pl	91
pt	97
ro	103
ru	109
sv	115
tr	121

Consumer information

Please read carefully

Patient Information Leaflet

- Read this leaflet carefully because it contains important information for you.
- Keep this leaflet. You may need to read it again.
- Ask your dentist if you need more information or advice.
- You must contact your dentist if any of the side effect describe further below worsen or do not improve.

This Patient Information Leaflet includes important information for you to consider and discuss with your healthcare provider prior to treatment with creos™ mucogain, including the benefits and risks of the treatment and any precautions which should be taken.

Soft Tissue Matrix (STM)

creos™ mucogain

What is creos™ mucogain?

Your healthcare provider determined a need for a soft tissue regeneration (filling and/or adding volume in the mouth around teeth and/or implants) and for stabilization of oral wounds.

creos™ mucogain is a collagen resorbable matrix. It is composed of purified natural protein collagen sourced from porcine tissues (pigs).

creos™ mucogain has sponge-like structure that intermingles with new (soft) tissue and is then included in the natural remodeling process of soft tissue.

Why is a Soft Tissue Matrix needed?

creos™ mucogain provides a supporting structure, matrix, for new soft tissue formation around teeth and/or implants and supports the damaged and/or missing tissue.

In this way, health care provider does not need to cut out your own tissue (so called autologous graft provided from patients own body) but uses creos™ mucogain instead.

Does the matrix have to be removed in a second surgery?

No. Over an 8-week period time, the collagen matrix is widely degraded by the body. As the new tissue is forming, the collagen matrix disappears. The matrix is completely absorbed by the body. A second operation to remove the matrix is therefore not needed.

creos™ mucogain consists of tissue of animal origin. How safe and biocompatible is it?

Starting materials for the production of creos™ mucogain are harvested in EU certified facilities under strict veterinary controls from veterinary certified animals that are declared fit for human consumption. The biocompatibility of creos™ mucogain has been demonstrated through testing and is due to the utilization of effective purification and processing methods during manufacturing. It complies with all relevant regulations ensuring its safety and effectiveness.

Composition of the device

creos™ mucogain consists of: collagen and elastin fibers (84.7 %). In addition, it contains

13 % moisture (H₂O), 2.3 % sodium chloride (NaCl), trace amounts of hydrochloric acid (HCl) and trace amounts of ethanol (EtOH).

When shall creos™ mucogain not be used?

The following are contraindications for creos™ mucogain:

- creos™ mucogain is not indicated for exposed healing.
- creos™ mucogain should not be used in infected areas. An active infection of the implant site may lead to an increased biodegradation of creos™ mucogain which may lead to an early loss of its matrix function. In such cases, the infection or inflammation should be adequately treated first.
- Adverse reactions during the use of collagen based matrix during dental soft tissue regeneration procedures have not been observed. However, allergic reactions towards porcine collagen may not be totally excluded. The use of creos™ mucogain in patients with known sensitivity to porcine-derived materials and collagen is not recommended.
- creos™ mucogain has not been evaluated in pediatric/adolescent patients and is not recommended for use in children. Routine treatment is not recommended until the end

of the juvenile jaw bone growth phase has been properly documented.

- Clinical data on the use of creos™ mucogain during pregnancy or lactation is currently not available. As a safety precaution, it is not recommended to use creos™ mucogain in women who are pregnant or breast feeding.
- creos™ mucogain should not be used to increase keratinized tissue (tissue that holds the teeth and implants in place and protects the roots) unless a minimal amount of this tissue is already present and preserved in the recipient area before surgical preparation and treatment.

Where are precautions necessary?

In patients with constrained healing capacity (e.g. due to metabolic disorders, cancer treatments, heavy smoking) creos™ mucogain should be used with special caution and closer patient monitoring during postoperative care.

How do I attend my wound?

A bacterial inflammation of the gum and jawbone is one of the rare complications that can occur after creos™ mucogain treatment. The good news is that by maintaining good oral hygiene, you can generally prevent it by yourself!

Providing good preventive care is, therefore, the best thing you can do to take care of your wound.

In the first few days after surgery, you must protect the sensitive wound, so your dentist may have recommended to not brush the operated area and/or area with the new implant and may have instead suggested a rinsing regime which may also include an antiseptic mouthwash such as containing chlorhexidine.

To further encourage good healing during the first few days you should also avoid spicy foods as well as coffee, black tea, and alcohol, and allow your food to cool prior to eating. Furthermore, smoking is not recommended during the first days after surgery.

Depending on the surgical protocol that was chosen by your dentist to be best suitable for you, eating and brushing recommendations may vary.

Are there any undesirable side effects possible?

The possible side effects related to the surgical procedure for creos™ mucogain can be: infection, wound swelling, bleeding, local inflammation, matrix exposure and/or pain.

Please seek medical advice if you experience any of the forementioned symptoms.

Magnetic Resonance Imaging (MRI) with creos™ mucogain

Since creos™ mucogain does not contain any metallic components, there is no impairment of creos™ mucogain during an MRI scan.

For a patient in the European Union and in countries with an identical regulatory regime:

If, during the use of this device or as a result of its use, a serious incident has occurred in relation to creos™ mucogain, please report it to the manufacturer and to your national authority. The contact information for the manufacturer of this device to report a serious incident is as follows:

Phone: +49 2407 56440

Email: vigilance@matricel.com

According to the European Medical Device Regulation (MDR (EU) 2017/745) a Summary of Safety and Clinical Performance document (SSCP) is available for the creos™ mucogain product. The SSCP can be obtained at the following website:
<https://ec.europa.eu/tools/eudamed>¹

¹Website available upon launch of the European Database on Medical Devices (EUDAMED).

In the meantime, the SSCP can be made available on request.

This information sheet is available for download at: <http://www.matricel.com/patient>

Implant card

Together with this Patient Information Leaflet you have received a Patient Implant Card (PIC) (or will receive from your healthcare provider soon after the treatment).

This card, similar to a business card, contains all relevant information related to your individual creos™ mucogain.

You can keep this card in your wallet to always have that information with you.

Symbols used on the Implant Card

Symbol	Description
--------	-------------

Patient Name or patient ID

(This is filled in by the implanting healthcare institution, either printed, handwritten, by placing a label e.g.)

Date of implantation

(This is filled in by the implanting healthcare institution, either printed, handwritten, by placing a label e.g.)

Name and Address of the implanting healthcare institution/provider

(This is filled in by the implanting healthcare institution, either printed, stamped, handwritten, by placing a label e.g.)

Name and Address of the manufacturer of the soft tissue matrix.

A website of the manufacturer where patients can get more information about the implant.

Device Name

This symbol is used to indicate the device name.

Lot Number

This is the manufacturer's unique lot number of your implant.

Use-by Date

This is the date after which the medical device may no longer be implanted. Indicated on the Documentation label that is placed on the Implant card.

Unique device identification (UDI) in automatic identification and data capture (AIDC) format.

Manufacturer's information coded in here in a digitally readable format.

A human readable format of the UDI information.

Contains biological material of animal origin.

This leaflet is Rev 0623 and was issued on 2023-06-29.

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© Nobel Biocare Services AG, 2023 All rights reserved. Nobel Biocare, the Nobel Biocare logotype and all other trademarks used in this document are, if nothing else is stated or is evident from the context in a certain case, trademarks of Nobel Biocare.
Matricel GmbH
Nobel Biocare AB